

Parameter

Tip Hardness	30HRC
Amplitude	20-60µm
Treatment Depth	14mm
Power Input	100V-240V 50/60HzL.2A(Max)
Input Power	48VA
Tip Offset	1µm100µm Root Canal
Tip Offset	<1mm
Half Offset Force	0.1N2N
Operating Frequency	2040kHz
Circuit Board Fuse	F1.6AL250V
Main Unit Fuse	T0.5AL250V
Dimensions	290mmx168mmx155mm



Xpedent Scaler Machine

SHEL-P



DIFFUSION
LASER & TECHNOLOGY
www.diffusion-laser.com

Résidence de Virgile - B4
18 rue de la Duchesse de Chartres
60 500 Vineuil-Saint-Firmin
Mobile: +33 (0) 611 417010
mail: contact@nf-diffusion.com



DIFFUSION
LASER & TECHNOLOGY
www.diffusion-laser.com

SHEL-P



Autoclave box × 1




Pièce a Main X2



Clé pour insertx 1



Pedal × 1



eAx2 ePx1 ePSx1 ePL3x1 eICP SETx1

6 INSERTS



Connection Externe de L'eau

Réservoir d'eau ou de solution

Control du débit du Spray

Sortie du cordon

Freely switched built-in and external water supply

Ecran control Tactil et Nettoyable Etanche

Eclairage Led

Branchement de la Pièce a Main

